


BACTRACK® S50  
SELECT

# ALCOHOL DETECTOR

TESTS BREATH FOR THE PRESENCE  
OF ALCOHOL AND ESTIMATES BAC



- Simple One-button Operation
- Get Immediate Results
- Comes with Hygienic Mouthpieces



50 Mouthpieces Included



Uses Advanced  
FlowCheck™ Sensing  
Technology

**Great for personal use, schools,  
businesses, clinics, and more!**

**0.00-0.40%**  
BAC TESTING RANGE

Meets Department of Transportation  
Specifications for Breath Alcohol Screeners

SELECT SERIES

**BACTRACK® S50**  
SELECT

# ALCOHOL DETECTOR

TESTS BREATH FOR THE PRESENCE OF ALCOHOL AND ESTIMATES BAC

The BACTRACK® Select S50 quickly and easily screens for the presence of alcohol. Simply blow into the mouthpiece and a BAC estimate is displayed within seconds.



## Easy-To-Use Hygienic Mouthpieces

Great for multiple users.

## Large LED Display

Easily viewable in daylight or at night.

### PRODUCT SPECIFICATIONS

Detection Range .....0.00 – 0.40 %BAC  
Power Supply .....Two AA batteries, included  
Warm Up Time .....20-30 seconds  
Response Time .....5 seconds  
Operating Temperature ....50-104 °F (10-40 °C)

### KIT INCLUDES

- BACTRACK® Select S50
- One-year limited warranty
- Instruction manual
- 2 AA batteries
- 50 Mouthpieces
- Soft carrying pouch

### INTENDED USE

This device is intended to measure alcohol in human breath. Measurements obtained by this device are used in the diagnosis of alcohol intoxication. Estimates and results may vary. This product cannot determine the subject's exact blood alcohol concentration (BAC) and should not be used for such purposes. Read the enclosed instructions first to ensure proper use, accuracy limitations and calibration requirements. It is very important to have the instrument calibrated as recommended in the instructions. Failure to do so may result in inaccurate readings. Do not rely upon results to determine level of intoxication or your ability to safely operate motor vehicles or heavy machinery.

[WWW.BACTRACK.COM](http://WWW.BACTRACK.COM)

© 2009 KHN Solutions LLC  
300 Broadway, Suite 26  
San Francisco, CA 94133  
877.334.6876

Designed in California  
Assembled in Korea

